

**Product**  
**Continuous oven**  
**Type**  
**VDU 55/40/200-100°C**

**Application**

Drying of pure-water impregnated filtration cartridge

**Branch**

Pharmaceutical industry

**Technical Data**

Nominal temperature: 100°C  
Heating capacity: 18 kW  
Supply voltage:  
3/N/PE AC 400 V ± 10%, 50/60 Hz

**Inner dimensions:**

Width 1000 mm  
Height 1500 mm  
Depth 1000 mm

**Outer dimensions**

Width 1600 mm  
Height 3300 mm  
Depth 1340 mm

**System description**

- Outer housing made of stainless steel (material no. 1.4301), polished surface (grain 400)
- Inner housing made of stainless steel (material no. 1.4301) in compliance with GMP
- Charging via trolley and auxiliary trolley
- SIMCON control
- Hot water tempering device 140°C

**Special features**

- Surface roughness < 0.8 µm
- Installation in a clean room of Class ISO 7 (DIN EN ISO 14644-1)

